

Baxter

The new modality “Expanded Hemodialysis,
HDx” for the Home dialysis patients to clear
Large Middle Molecules

Frans Richter
Medical Manager Renal Care
Baxter Benelux



Theranova: A New Dimension


By the end of this session you will

- Understand the evolution of this dialyzer
- Understand how Theranova and HDx are the next step Hemodialysis

Agenda


- Renal Transport Middle Molecules
- HDx therapy with the new MCO Theranova membrane
- Summary

We have not solved cardiovascular disease



Adapted from Sarnak MJ et al. Hypertension 2003;42:1050-1065

Survival in Dialysis?




Large Middle Molecules:

- Meet the criteria for an uremic toxin
- Are middle-molecules which currently have limited removal by standard low and high flux HD

Molecular weights >15kDa

Renal transport of middle molecules

Dependence on solute size and molecular weight





Weekly solute clearance as a function of molecular weight

Cheung AK. *Blood Purif* 1994;12:42-53

Uremic retention solutes of the Middle Molecule group indicating solutes ≥ 15 kDa and ≥ 20 kDa in size

α_1-Acid glycoprotein	Desacylghrelin	Myoglobin
α_1-Microglobulin	Dinucleoside polyphosphates	Neuropeptide Y
Adiponectin	β -Endorphin	Orexin A
Adrenomedullin	Endothelin-1	Orexin B
AGEs	Fibroblast growth factor-23	Osteocalcin
Angiogenin	κ-Ig light chain	PTH
Angiotensin A	λ-Ig light chain	Pentraxin-3
AOPPs	Ghrelin	Prolactin
Atrial natriuretic peptide	Glutathion, oxidized	Resistin
Basic fibroblast growth factor	Guanylin	Retinol binding protein
β_2 -Microglobulin	Hyaluronic acid	Soluble intracellular adhesion molecule-1
β-Trace protein	IGF-1	Soluble TNF receptor 1
Calcitonin	IL-1β	Soluble TNF receptor 2
Calcitonin-gene-related peptide	IL-6	Substance P
Cholecystokinin	IL-8	TNF-α
Clara cell protein (CC16)	IL-10	Uroguanylin
Complement factor D	IL-18	Vascular endothelial growth factor
Cystatin C	Leptin	Vasoactive intestinal peptide
Degranulation inhibiting protein I	Methionin-enkephalin	Visfatin
Delta-sleep inducing peptide	Motilin	YKL-40

 Molecular weight 15-20 kDa

23  Molecular weight 20 kDa or greater

Vanholder et al. *Kidney Int* 2003;63:1934–43
 Durantou et al. *J Am Soc Nephrol* 2012;23:1258–70
 Neiryneck et al. *Int Urol Nephrol* 2013;45:139–50
 Chmielewski et al. *Sem Nephrol* 2014;34:118–34
 Hutchison C. ERA-EDTA 2016 Presentation

Large Middle Molecule uremic solutes and potential impact on health


α_1 -Acid glycoprotein	[43 kDa]	Inflammatory processes (2)
α_1 -Microglobulin	[33]	Restless leg syndrome (6)
Adiponectin	[30]	Unclear, low concentration is likely a reflection of abdominal fat/nutritional status (1)
AGEs	[variable]	Inflammatory processes, atherosclerosis progression (1)
Basic fibroblast growth factor	[24]	
β -trace protein	[26]	
Complement factor D	[24]	Chronic inflammation, atherosclerosis progression (1)
Fibroblast growth factor-23 (FGF-23)	[32]	Uremic cardiomyopathy (4)
κ -Ig light chain	[22.5]	Impaired immune defense, nephrotoxicity (1)
λ -Ig light chain	[45]	Impaired immune defense, nephrotoxicity (1)
Hyaluronic acid	[variable]	
Interleukin-1 β (IL-1 β)	[32]	Inflammation, cognitive dysfunction (5)
Interleukin-6 (IL-6)	[24.5]	Inflammation, atherosclerosis progression, anemia (1); cognitive dysfunction (5)
Interleukin-18 (IL-18)	[20]	
Pentraxin-3	[40]	Endothelial dysfunction (1)
Prolactin	[22]	Sexual abnormalities, anemia, endothelial dysfunction, arterial stiffness (1)
Retinol binding protein	[21]	Oxidative stress, atherosclerosis progression (1)
Soluble TNF receptor 1	[30]	Prolong half-life of TNF α and Increase its cytotoxicity
Soluble TNF receptor 2	[40]	Prolong half-life of TNF α and Increase its cytotoxicity
TNF- α	[26]	Coagulation disorders, insulin resistance, endothelial dysf. wasting (1); cognitive dysf. (5)
Vascular endothelial growth factor	[34]	
Visfatin	[55]	Endothelial damage, inflammation, plaque destabilization (1)
YKL-40	[40]	Associated with standard inflammatory parameters (3)

(1) Chmielewski et al. *Sem Nephrol* 2014;34:118–34; (2) Lisowska Myjak B. *Nephron Clin Pract* 2014;128:303-11; (3) Okyay GU et al. *Ther Apher Dial.* 2013;17:193-201; (4) Grabner A et al. *Curr Opin Nephrol Hypertens* 2016;25:314-24; (5) Watanabe et al *NeuroToxicology* 2014;44:184-93; (6) Sakurai K. *Blood Purif* 2013;35(Suppl1.):64-8;


Agenda

- Renal Transport Middle Molecules
- **HDx therapy with the new MCO Theranova membrane**
- Summary

Current therapies have limitations in the clearance of large uremic solutes (large middle molecules)





Expanded hemodialysis, HDx: A step closer to the native kidney in which molecules are cleared, simpler than HDF



Membrane pore size distribution is determined by the spinning process conditions


Innovation on the spinning process and speed of porous hollow fibers




Phase Inversion Process

Kim HJ. *Jpn J Appl Phys* 2016;55 06GH06

non-uniform pore sizes (highflux)





pore size distribution TheraNova




MCO membrane characteristics support internal convection

Inner fiber diameter



Polyflux



Theranova MCO

- Increased performance as sum of:
- Increase of membrane permeability
 - Reduction of inner diameter, more fibers per dialyzer
- Optimum (Performance, Process Stability): 180 / 35 μm



The Rise of Expanded Hemodialysis

Claudio Ronco^{a, b}

^aDepartment of Nephrology Dialysis and Transplantation, St. Bortolo Hospital, and ^bInternational Renal Research Institute, Vicenza, Italy

Expanded Hemodialysis

Innovative Clinical Approach in Dialysis
Editor(s): Ronco C. (Vicenza)
Status: available
Publication year: 2017



Membrane as barrier to microbiological contaminants in dialysis fluid – Simulating clinical conditions, a study by the Ghent team

LPS detection by LAL test (EU/ml)

mean ± standard deviation [range] {number of repeats out of 6 above detection limit}

Product	Dialysate	Blood side	
	(natural ET challenge)	at start	after 1 hour
Low-flux (Polyflux L)	8.6 ± 5.8 [3.5-19.1]	0.004 ± 0.000 [0.004-0.004] {0}	0.005 ± 0.002 [0.004-0.008] {1}
High-flux (Revaclear)	12.2 ± 12.2 [3.6-33.7]	0.005 ± 0.002 [0.004-0.008] {2}	0.005 ± 0.001 [0.004-0.008] {2}
Medium cut-off (Theranova)	8.3 ± 2.4 [6.0-11.8]	0.004 ± 0.000 [0.004-0.004] {0}	0.006 ± 0.004 [0.004-0.014] {3}
High cut-off (Theralite)	8.9 ± 7.4 [3.2-22.5]	0.004 ± 0.001 [0.004-0.006] {1}	0.007 ± 0.005 [0.004-0.016] {5}

bacterial endotoxin (*lipopolysaccharide, LPS*)

detection limit: 0.005 EU/ml

Endotoxin retention for dialysis membranes of different permeability profile


Membrane type ¹	Pore radius [nm] (mean ± SD)	Endotoxin LRV ² (mean ± SD)
Low-flux	3.1 ± 0.2	2.8 ± 0.2
High-flux	4.5 ± 0.2	3.3 ± 0.3
MCO 4	6.5 ± 0.2	3.5 ± 0.1
HCO	10 ± 2	3.3 ± 0.5

¹ All membranes made of same material, having similar wall structure

² Tested with LPS from *E.coli* O55:B5 (N=3); LRV = Logarithmic Retention Value

Ref: Hulko et al. ERA-EDTA 2015 abstract FP516

Reported albumin removal in HDF and HDx treatments



HDF data are obtained in different studies using a variety of high-flux dialyzers, different dilution modes (post-, pre-, mid-), and different convective flow rates


Theranova Dialyzer in HD mode

N = 19

$Q_B = 301 \pm 22$ ml/min

T = 4.0 h

Pre- to post-dialysis reduction in plasma level



High-flux HD = HD by FX Cordiax 80 dialyzer

Bars indicate mean and SD

Post-dialysis data are corrected for hemoconcentration

Kirsch et al. *Nephrol Dial Transplant* 2017;32:165-72

HDx performance versus HDFoL and high-flux HD


N = 19

$Q_B = 400$ ml/min

T = 4.4 ± 0.3 h

HDF $V_{CONV} = 24$ L

Pre- to post-dialysis reduction in plasma level



hvHDF = high volume HDF by FX Cordiax 800; high-flux HD = HD by FX Cordiax 80

Bars indicate mean and SD (except for alpha1microglobulin data that are presented as medians)

Post-dialysis data are corrected for hemoconcentration

Statistics by a mixed model with fixed effects of period and study dialyzer type, and the random effect of subject.

What could improve with increased removal of large MM?

- Chronic inflammation
- Cardiovascular disease
- Infections (secondary immunodeficiency)

- Quality of life

Agenda

- Renal Transport Middle Molecules
- HDx therapy with the new MCO Theranova membrane
- **Summary**

Published clinical studies of MCO Membrane to date

Author	Study	Patients (n)	Abstract	Paper
Kirsch AH	Clearance	39		Neph Dial Trans (2017)
Zickler D	MCO Ci	48		PLOS One 2017
Belmouaz M	HDF to HDx	10		Clin Nephrology 2017
Teatini U	HDF vs HDx	8	ERA-EDTA; 2017	
Schepers E	<i>(In vitro)</i>		ERA/EDTA 2017	
Koball S		34	ASN 2017	
Cantaluppi	HF vs HDx	14	EDTA 2018	
Celik L	HF vs HDx	8	EDTA 2018	
Giuseppe G	HF to HDx	14	EDTA 2018	
Santivañez J	HDF vs HDx	8	EDTA 2018	

HDx/MCO Membrane Studies

Study	Key endpoints	Comparator	Study duration per subject	# Subjects
REMOVAL-HD PI: Hutchinson, ANZ multi-center [ACTRN12616000804482]	Changes in S-albumin and FLC levels over 6 months; Several exploratory measures	Standard of care high-flux HD	8 months	85
PI: Koball, Rostock, Germany [DRKS00011638]	Albumin-binding capacity	HDF	4 weeks	30
REMOG PI: Rosenkranz, Graz, Austria [DRKS00012359]	Change in FLCs and S-albumin, various MMs and cytokines, calcification propensity, patient-reported outcomes	HDF	6 months	30
PI: Bridoux, Poitiers, France [NCT03211676]	MMs and inflammatory markers, several secondary measures	HD	6 months	40
ModuVas PI: Zickler, Berlin, Germany [NCT03104166]	Vascular calcification in human cell culture model, PWV and T50 calcification propensity test	HD	7 months	48
PI: Cozzolino, Milan, Italy [NCT03169400]	Vascular calcification and oxidative stress in rat cell culture model	HD	6 months	20
PI: Juillard, Lyon, France [NCT03137056]	I: Proteome of depuration, inflammatory markers II: Effect on albumin isoforms	I: HD II: HD/HDF	I: 16 weeks II: 24 weeks	I: 20 II: 20
MCO-IF PI: Schmaderer, Munich Germany [NCT03270371]	Change in inflammatory score ; calcification biomarkers, calcification propensity	HD	8 months	50
PI: Caldin da Silva, Sao Paolo, Brazil [NCT03274518]	Medium molecule removal; intradialytic hemodynamics; fluid status	HDF	2½ months	16

Key Takeaways

HDx enabled by *Theranova is*

- 1) The next step in dialyzer evolution
- 2) The new MCO membrane with larger pores and increased selectivity that provides in HD mode greater removal of large middle molecules than conventional high-flux membranes
- 3) A pure HD dialyzer!

Please keep
in mind:

